# SANITAS

## SBF 14

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## ENGLISH

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

# ! 1. Important notes – keep for later use!

### 1.1 Safety notes:

- **The scales must not be used by persons with medical implants (e.g. heart pacemakers). Otherwise their function could be impaired.**
- Do not use during pregnancy.
- Caution! Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- Batteries are highly dangerous if swallowed. Keep batteries and scale out of reach of small children. If batteries are swallowed, get medical help immediately.
- Keep children away from packing materials (danger of suffocation).
- Batteries must not be charged or reactivated by other means, taken apart, thrown into a fire or short-circuited.

### 1.2 General notes:

- Not intended for medical or commercial use.
- Please note that measuring tolerances are possible due to technical reasons because the scale has not been calibrated.
- Ages 10 ... 100 years and height settings from 100 ... 220 cm (3'-03"–7'-03") can be preset. Load capacity: max. 180 kg (396 lb, 28 st.). Result in steps of 100g (0.2 lb, 0.01st). Measurements of body fat, body water and muscle density in steps of 0.1%.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm". To the rear of the scales, there is a switch where you can select 'pound' and 'stone' (lb, st).
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- The scale should be cleaned occasionally with a damp cloth. Do not use abrasive detergents and never immerse the scale in water.
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- Repairs may only be carried out by customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.

## 2. Operation

### 2.1 Inserting batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity. If the scale fails to operate, remove the battery completely and insert it again. Your scale is equipped with a "replace battery" indicator. When the scale is operated with a battery which is too weak, "L<sub>B</sub>" appears in the display field and the scale automatically switches off. In this case the battery must be replaced (2 x 3V lithium battery CR2032).

### 2.2 Setting user data

"SET"

User data	Set values
Memory position P	1 to 10
Height cm	100 to 220 cm (3'-03" to 7'-03")
Age	10-100
Sex	male  , female 
Degree of activity R	1-5

Use the  and  buttons to select and confirm with "SET".

## Degrees of activity

Degrees of activity	Physical activity
1	None.
2	Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises).
3	Medium: Physical effort for 30 minutes at least 2 to 4 times a week.
4	High: Physical effort for 30 minutes at least 4 to 6 times a week.
5	Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily.

## 3. Operation

### 3.1 Measuring weight only

- Stand on the scale without moving about and distribute your weight evenly on both legs. The scale immediately begins to measure your weight. Your measured weight is displayed shortly after this (Fig 1).
- The scale switches off a few seconds after you step down from it.

Fig. 1

### 3.2 Measurement of weight, body fat percentage and body water content, muscle percentage and bone mass

- Switch on the scale (tap the standing surface).
- Select the memory preset where your basic personal data are stored by pressing the  $\checkmark$  or  $\wedge$  button repeatedly. The data are displayed consecutively until "0.0" appears.
- Step onto the scale barefoot and make sure you're standing on both electrodes. First your body weight is determined and displayed (Fig. 1).
- Remain on the scale and avoid moving; your body fat is being analysed. This may take a few seconds.

#### The following data are displayed:

- Weight in kg (lb, st) (Fig. 1)
- Body fat percentage **FAT** in % (Fig. 2)
- Water percentage **BW** in % (Fig. 3)
- Muscle percentage **MUS** in % (Fig. 4)
- Bone mass in **BON** in kg (Fig. 5)
- Now all measured values are displayed consecutively and the scale switches off.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

### 3.3 Tips for using the scale

#### Important when measuring body fat/body water/muscle percentage:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Your feet, legs, calves and thighs must not touch each other, as otherwise the measurement cannot be performed correctly.
- Stand still during the measurement.
- Wait several hours (6-8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

#### The measurement is not reliable for:

- Children under approx. 10 years of age.
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis.
- Persons taking cardiovascular medication.
- Persons taking vasodilating or vasoconstricting medications.
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

## 4. Evaluation of results

### Body fat percentage

The following body fat percentages are guide values (contact your physician for further information).

#### Man

Age	low	normal	high	very high
10-14	<11	11-16	16,1-21	>21,1
15-19	<12	12-17	17,1-22	>22,1
20-29	<13	13-18	18,1-23	>23,1
30-39	<14	14-19	19,1-24	>24,1
40-49	<15	15-20	20,1-25	>25,1
50-59	<16	16-21	21,1-26	>26,1
60-69	<17	17-22	22,1-27	>27,1
70-100	<18	18-23	23,1-28	>28,1

#### Woman

Age	low	normal	high	very high
10-14	<16	16-21	21,1-26	>26,1
15-19	<17	17-22	22,1-27	>27,1
20-29	<18	18-23	23,1-28	>28,1
30-39	<19	19-24	24,1-29	>29,1
40-49	<20	20-25	25,1-30	>30,1
50-59	<21	21-26	26,1-31	>31,1
60-69	<22	22-27	27,1-32	>32,1
70-100	<23	23-28	28,1-33	>33,1

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

### Body water percentage

The body water percentage is normally within the following ranges:

#### Man

Age	poor	good	very good
10-100	<50	50-65	>65

#### Woman

Age	low	normal	high
10-100	<65	45-60	>60

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

### Muscle percentage

The muscle percentage % is normally within the following ranges:

#### Man

Age	low	normal	high
10-14	<44	44-57	>57
15-19	<43	43-56	>56
20-29	<42	42-54	>54
30-39	<41	41-52	>52
40-49	<40	40-50	>50
50-59	<39	39-48	>48
60-69	<38	38-47	>47
70-100	<37	37-46	>46

#### Woman

Age	low	normal	high
10-14	<36	36-43	>43
15-19	<35	35-41	>41
20-29	<34	34-39	>39
30-39	<33	33-38	>38
40-49	<31	31-36	>36
50-59	<29	29-34	>34
60-69	<28	28-33	>33
70-100	<27	27-32	>32

### Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. The bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Please note that this scale does not identify the calcium content of the bones; instead, it calculates the weight of all the components that make up the bones (organic matter, inorganic matter and water). It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

#### IMPORTANT:

Please do not confuse bone mass with bone density.

The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

## 5. Incorrect measurement

If the scale detects an error during weighing, "FFFF" or "Err." appears in the display.

If you step onto the scale before "0.0" appears in the display, the scale will not operate properly.

Possible causes of errors:	Remedy:
- The maximum load-bearing capacity of 180 kg (396 lbs) was exceeded.	- Only weigh the maximum permissible weight.
- The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	- Repeat weighing barefoot. - Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
- Your body fat lies outside the measurable range (less than 5% or greater than 55%).	- Repeat weighing barefoot. - Slightly moisten the soles of your feet if necessary.
- Scale is on carpeting.	- Place the scale on a level and solid surface.

## 6. Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

### Note

The codes below are printed on batteries containing harmful substances: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life.

Please dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please contact the appropriate local authorities.

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